



MENU- Full List of Settings

- O IRIS ELC / ALC
- O FOCUS ADJ OFF / ON
- O EXPOSURE
 - BRIGHTNESS 12|-----|
 - SHUTTER AUTO _|
 - MODE NORMAL/DEBLUR
 - RETURN _|
 - SENS UP OFF / X2 / X4 / X8 ... X64
 - AGC 2|-----|
 - RETURN _|
- O BACKLIGHT OFF / HLC _| / BLC _| / WDR _|
 - HLC
 - LEVEL 10|-----|
 - COLOR BLK / WHT / YEL / CYN / GRN / MAG / RED / BLU
 - RETURN _|
 - BLC
 - H-POS 7 0-20
 - V-POS 6 0-20
 - H-SIZE 5 0-20
 - V-SIZE 5 0-20
 - RETURN _|
 - WDR
 - WEIGHT MIDDLE / HIGH / LOW
 - RETURN _|
- O DAY&NIGHT COLOR / B&W _| / EXTERN _| / AUTO _|
 - B&W
 - IR LED OFF / ON _|
 - ANTI-SAT. 10|-----|
 - RETURN _|
 - EXTERN
 - IR LED OFF / ON _|
 - ANTI-SAT. 10|-----|
 - EXTERN SW HIGH / LOW / MIDDLE
 - D > N THRES 13|-----|
 - N > D THRES 7|-----|
 - DELAY LOW / MIDDLE / HIGH
 - RETURN _|
 - AUTO
 - IR LED OFF / ON _|
 - ANTI-SAT. 10|-----|
 - AGC THRES 10|-----|
 - AGC MARGIN 10|-----|
 - DELAY LOW / MIDDLE / HIGH
 - RETURN _|
- O COLOR
 - AWB PRESET v / MANUAL _| / AUTO / AUTOext
 - MANUAL
 - COLORTEMP 3000K / 5000K / 8000K
 - R-GAIN 10|-----|
 - B-GAIN 10|-----|
 - COLOR GAIN 2|-----|
 - RETURN _|
- O DNR MIDDLE
- O IMAGE
 - SHARPNESS 5|-----|
 - GAMMA 0.55 / 0.65 / 0.75 / 0.45
 - MIRROR OFF / ON
 - FLIP ON/OFF
 - D-ZOOM 1.0x -> 16.0x: 0.1x
 - ACE OFF / LOW / MIDDLE / HIGH
 - DEFOG OFF / ON _|
 - MODE MANUAL / AUTO
 - WEIGHT LOW / MIDDLE / HIGH
 - RETURN _|
 - SHADING
 - WEIGHT 1% - 100%
 - RETURN _|
 - PRIVACY
 - BOX OFF / ON _|
 - ZONE NUM 0-15
 - ZONE DISP OFF / ON
 - H-POS 0-60



Pharos HD

Menu Guide

Issued by: Elliott Lasher

Eff. Date:

Rev: 0

Pg. 2 of 4

```

V-POS          0-34
H-SIZE         0-60
V-SIZE         0-34
Y LEVEL 10 |----|----|
CB LEVEL 10 |----|----|
CR LEVEL 10 |----|----|
TRANS          3 / 2 / 1 / 0
RETURN        _|
POLYGON OFF / ON _|
ZONE NUM      0-15
ZONE DISP     OFF / ON
POS1-X        0-120
POS1-Y        0-68
POS2-X        0-120
POS2-Y        0-68
POS3-X        0-120
POS3-Y        0-68
POS4-X        0-120
POS4-Y        0-68

Y LEVEL 10 |----|----|
CB LEVEL 10 |----|----|
CR LEVEL 10 |----|----|
TRANS          3 / 2 / 1 / 0
RETURN        _|
O  DIS        RETURN        _|
                OFF / ON _|
                RANGE        10% / 20% / 30%
                FILTER        MIDDLE / HIGH / LOW
O  iMOTION    AUTO C        FULL / OFF / HALF
                OFF / ON _|
                DET WINDOW    _|
                WINDOW ZONE0 / 1 / 2 / 3
                WINDOW USE    ON / OFF
                DET H-POS      1-60
                DET V-POS      1-34
                DET H-SIZE     1-60
                DET V-SIZE     1-34
                RETURN        _|
                SENSITIVITY    5 |----|----|
                ESTIMATE OFF / ON
                QUCIK ZOOM     OFF / ON _|
                ZOOM SPEED     MIDDLE / HIGH / LOW
                TRACKING       OFF / ON
                REPEAT         OFF / ON
                RETURN        _|
                MOTION OSD     OFF / ON
                TEXT ALARM     OFF / ON
                SIGNAL OUT     OFF / ON
                RETURN        _|
O  SYSTEM     _|
                COM           _|
                CAMID          1-255
                BAUDRATE       115200
                RETURN        _|
                IMAGERANGE     FULL / COMP / USER _|
                OFFSET         3 |----|----|
                RETURN        _|
                COLORSPACE     HD-CbCr / YUV / SD-CbCr
                FRAMERATE      1080-60p / 720-30p / 720-60p / 1080-60i / 1080-30p
                FREQ           60Hz / 50Hz
                CVBS           OFF / ON
                COLOR BAR      OFF / ON
                LANGUAGE       ENG / CHN / CHN(S) / JPN / KOR
                CAM TITLE      OFF / ON _|
                                0000000
                                -----
                                U, D to change value
                                R, L to change location
                                Return to enter
                RESET          ONv
O  EXIT       RETURN        _|

```



Pharos HD

Menu Guide

Issued by: Elliott Lasher

Eff. Date:

Rev: 0

Pg. 3 of 4

Refer to the manuals and the Operating Card for normal Pharos HD™ use.

Access the menu by pressing and holding the ON/OFF button, followed by pressing the center White Balance button. Access will be apparent by the strobing of a red/green LED by the ON/OFF button. Any connected display will have a list of settings in white letters. Use the Zoom and Brightness buttons as arrows and the White Balance button as an enter button. The current selected option will be shown with yellow letters.

Default Menu List and important setting options

IRIS ELC

The iris options are set between an electronic iris, and an automatic iris. Since the Pharos HD uses a manual lens with a fixed iris, the ALC compensation is not recommended.

FOCUS ADJ OFF

The focus assist is not recommended. The image quality of even the unfocused image is enough to make the fine adjustments unclear.

EXPOSURE _|

The normal operation brightness buttons affect the first option in this menu. We recommend leaving the shutter, sensitivity and AGC settings. Changes will cause large reactions to the addition other surgical lights in and back out of the camera image.

Default brightness is set to 12, and default AGC is set to 2.

BACKLIGHT OFF

In the event of working with a known bright object on a specific part of the image, the backlight settings can be used to prevent the camera from changing brightness settings.

DAY&NIGHT COLOR

It is recommended to leave the camera set to the default option, Color.

COLOR _|

To increase or decrease color saturation, adjust the color gain. White balance settings can also be adjusted here. In order to use the white balance button effectively without having to set it, return to the preset mode and allow it to finish “pushing.”

Default White Balance mode is set to Preset, and default color gain is set to 4.

DNR MIDDLE

The difference between high middle and low is minimal in bright conditions.

IMAGE _|



Sharpness and Gamma have a decent impact on the image, and can be adjust for specific procedures. But the remaining options are recommended to leave at their default settings.

Default Sharpness is set to 5, and default Gamma is set to .55.

DIS OFF

DIS, or digital image stabilization, can be useful. The amount of the image periphery cut off to handle the image processing can be set to 10, 20 or 30%.

iMOTION OFF

This can be used to track motion of surgical tools.

SYSTEM _|

The camera can change resolutions, framerate, or display the serial number.

Default settings.

EXIT

Exit the menu by selecting the exit option. The diode will continue to blink red and green until the power button is hit once more. Normal button use will resume when the LED is green. The settings adjusted in this manner are saved to the device. Refer to the information provided for the default settings.